




CEMENTED FEMORAL STEM STANDARD (125° NECK ANGLE)									
Size	A	B				C			
	Stem Length (mm)	Horizontal Offset (mm)				Neck Length (mm)			
		-3.5	STD	+3.5	+7.0	-3.5	STD	+3.5	+7.0
00	127	34.2	37.0	39.9	42.8	32.6	36.1	39.6	43.1
01	149	42.2	45.0	47.9	50.8	37.1	40.6	44.1	47.6
02	149	42.2	45.0	47.9	50.8	37.1	40.6	44.1	47.6
03	149	42.2	45.0	47.9	50.8	37.1	40.6	44.1	47.6
04	149	42.2	45.0	47.9	50.8	37.1	40.6	44.1	47.6





CEMENTED FEMORAL STEM NARROW (125° NECK ANGLE)									
Size	A	B				C			
	Stem Length (mm)	Horizontal Offset (mm)				Neck Length (mm)			
		-3.5	STD	+3.5	+7.0	-3.5	STD	+3.5	+7.0
01	149	35.2	38.0	40.9	43.8	34.7	38.2	41.7	45.2
02	149	35.2	38.0	40.9	43.8	35.1	38.6	42.1	45.6
03	149	35.2	38.0	40.9	43.8	35.1	38.6	42.1	45.6
04	149	35.2	38.0	40.9	43.8	35.1	38.6	42.1	45.6



Libertas[®] Hip System

Rx only  Maxx Orthopedics, Inc.
 Norristown, PA 19403, USA

AJW Consulting GmbH
 Königsallee 106
 40215 Düsseldorf, Germany
 Phone: +49(0) 211 3013 2232
 Email: info@ajwtech.com

Carefully read all instructions and be familiar with the surgical techniques prior to use.

Please see the package inserts for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the LIBERTAS[®] Total Hip System.

CAUTION: THESE DEVICES ARE RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

WARNING: THE LIBERTAS[®] CEMENTED FEMORAL STEM IS INTENDED FOR CEMENTED USE ONLY.

LIBERTAS[®] Hip is manufactured by Maxx Orthopedics, Inc. LIBERTAS, LIBERTAS HIP, and Pursue Life are Registered Trademarks of Maxx Orthopedics, Inc.

©2020 Maxx Orthopedics. All rights reserved. Updated May 2020

Maxx Orthopedics Brochure No. MXO-LH00013e-R03