




CEMENTLESS FEMORAL STEM STANDARD (125° NECK ANGLE)													
Size	A	B						C					
	Stem Length (mm)	Horizontal Offset (mm)						Neck Length (mm)					
		-4.0	-3.5	STD	+3.5	+4.0	+7.0	-4.0	-3.5	STD	+3.5	+4.0	+7.0
00	115	41.7	42.2	45.0	47.9	48.3	50.8	34.8	35.3	38.8	42.3	42.8	45.8
01	130	42.2	42.7	45.5	48.4	48.8	51.3	34.8	35.3	38.8	42.3	42.8	45.8
02	140	43.2	43.7	46.5	49.4	49.8	52.3	34.8	35.3	38.8	42.3	42.8	45.8
03	145	43.7	44.2	47.0	49.9	50.3	52.8	34.8	35.3	38.8	42.3	42.8	45.8
04	150	44.7	45.2	48.0	50.9	51.3	53.8	34.8	35.3	38.8	42.3	42.8	45.8
05	154	45.2	45.7	48.5	51.4	51.8	54.3	34.8	35.3	38.8	42.3	42.8	45.8
06	160	45.7	46.2	49.0	51.9	52.3	54.8	34.8	35.3	38.8	42.3	42.8	45.8
07	165	46.7	47.2	50.0	52.9	53.3	55.8	34.8	35.3	38.8	42.3	42.8	45.8
08	170	47.2	47.7	50.5	53.4	53.8	56.3	34.8	35.3	38.8	42.3	42.8	45.8
09	180	48.2	48.7	51.5	54.4	54.8	57.3	34.8	35.3	38.8	42.3	42.8	45.8
10	189	49.2	49.7	52.5	55.4	55.8	58.3	34.8	35.3	38.8	42.3	42.8	45.8




CEMENTLESS FEMORAL STEM STANDARD (135° NECK ANGLE)													
Size	A	B						C					
	Stem Length (mm)	Horizontal Offset (mm)						Neck Length (mm)					
		-4.0	-3.5	STD	+3.5	+4.0	+7.0	-4.0	-3.5	STD	+3.5	+4.0	+7.0
00	115	35.2	35.5	38.0	40.5	40.8	43.0	34.8	35.3	38.8	42.3	42.8	45.8
01	130	35.7	36.0	38.5	41.0	41.3	43.5	34.8	35.3	38.8	42.3	42.8	45.8
02	140	36.7	37.0	39.5	42.0	42.3	44.5	34.8	35.3	38.8	42.3	42.8	45.8
03	145	37.2	37.5	40.0	42.5	42.8	45.0	34.8	35.3	38.8	42.3	42.8	45.8
04	150	38.2	38.5	41.0	43.5	43.8	46.0	34.8	35.3	38.8	42.3	42.8	45.8
05	154	38.7	39.0	41.5	44.0	44.3	46.5	34.8	35.3	38.8	42.3	42.8	45.8
06	160	39.2	39.5	42.0	44.5	44.8	47.0	34.8	35.3	38.8	42.3	42.8	45.8
07	165	40.2	40.5	43.0	45.5	45.8	48.0	34.8	35.3	38.8	42.3	42.8	45.8
08	170	40.7	41.0	43.5	46.0	46.3	48.5	34.8	35.3	38.8	42.3	42.8	45.8
09	180	41.7	42.0	44.5	47.0	47.3	49.5	34.8	35.3	38.8	42.3	42.8	45.8
10	189	42.7	43.0	45.5	48.0	48.3	50.5	34.8	35.3	38.8	42.3	42.8	45.8





CEMENTLESS FEMORAL STEM LATERAL (135° NECK ANGLE)													
Size	A	B						C					
	Stem Length (mm)	Horizontal Offset (mm)						Neck Length (mm)					
		-4.0	-3.5	STD	+3.5	+4.0	+7.0	-4.0	-3.5	STD	+3.5	+4.0	+7.0
00	115	42.2	42.5	45.0	47.5	47.8	50.0	39.0	39.5	43.0	46.5	47.0	50.0
01	130	42.7	43.0	45.5	48.0	48.3	50.5	39.0	39.5	43.0	46.5	47.0	50.0
02	140	43.7	44.0	46.5	49.0	49.3	51.5	39.0	39.5	43.0	46.5	47.0	50.0
03	145	44.2	44.5	47.0	49.5	49.8	52.0	39.0	39.5	43.0	46.5	47.0	50.0
04	150	45.2	45.5	48.0	50.5	50.8	53.0	39.0	39.5	43.0	46.5	47.0	50.0
05	154	45.7	46.0	48.5	51.0	51.3	53.5	39.0	39.5	43.0	46.5	47.0	50.0
06	160	46.2	46.5	49.0	51.5	51.8	54.0	39.0	39.5	43.0	46.5	47.0	50.0
07	165	47.2	47.5	50.0	52.5	52.8	55.0	39.0	39.5	43.0	46.5	47.0	50.0
08	170	47.7	48.0	50.5	53.0	53.3	55.5	39.0	39.5	43.0	46.5	47.0	50.0
09	180	48.7	49.0	51.5	54.0	54.3	56.5	39.0	39.5	43.0	46.5	47.0	50.0
10	189	49.7	50.0	52.5	55.0	55.3	57.5	39.0	39.5	43.0	46.5	47.0	50.0



Libertas® Hip System

Rx only  Maxx Orthopedics, Inc.
Norristown, PA 19403, USA

AJW Consulting GmbH
Königsallee 106
40215 Düsseldorf, Germany
Phone: +49(0) 211 3013 2232
Email: info@ajwtech.com

Carefully read all instructions and be familiar with the surgical techniques prior to use.

Please see the package inserts for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the LIBERTAS® Total Hip System.

CAUTION: THESE DEVICES ARE RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

LIBERTAS® Hip is manufactured by Maxx Orthopedics, Inc. LIBERTAS, LIBERTAS HIP, and Pursue Life are Registered Trademarks of Maxx Orthopedics, Inc.

©2020 Maxx Orthopedics. All rights reserved. Updated May 2020

Maxx Orthopedics Brochure No. MXO-LH00012e-R02